




AIVIA GO FOR AUTOMATIC EXTERNAL DEFIBRILLATORS (AED)



AED PUBLIC ACCESS PROGRAMS

Sudden cardiac arrest (SCA) is the first leading cause of death. In order to increase survival rates from sudden cardiac arrest, defibrillators' accessibility in public places is a solution recommended by medical authorities. Every AED has a visual indicator showing its operational status and AIVIA solutions complement AEDs to ensure their operability at all times.



AIVIA GO - SOLUTION FOR AED MONITORING

Aivia GO provides remote monitoring of the AED. Incident and events are retrieved via a software program and a wireless internet connection (4G), ensuring that the defibrillator remains fully operational and ready for use.

Aivia GO provides a direct voice communication with emergency services and a GPS tracking system allowing the emergency operator to provide clear instructions and assistance.

STANDARD FEATURES

EMERGENCY PHONE CALLS

- Emergency calling button to activate hands free telephone communication through LTE (911, 112 or any preset telephone number)
- 4G communications
- SIM Card support integrated

GPS TRACKING

- GPS tracking and immediate SMS notification when defibrillator is moved

AED STATUS MONITORING

- 24/7 Monitoring of defibrillator status, presence and temperatures
- Battery operated, autonomy up to 24 months

CONNECTIVITY

- Internet cloud-based software platform connectivity via 4G
- Automated SMS and emails for incidents and alerts



AIVIANET

A SOFTWARE MONITORING APPLICATION

- A multilingual version is available, requiring a username and password for access to the server
- Incidents are displayed with the possibility of verifying and logging facts
- One or several AIVIA GO can be monitored at any one time (status display and presence of the AED, control system, alarm...)
- AIVIANET can be programmed to send SMS/ Email alerts to one or more people should an incident occur
- Record keeping of past incidents
- Data interpretation by way of colour codes
- Geolocation of any AIVIA cabinet connected to the network



ACCESSORIES

TRANSPORT CASE



FASTENING STRAP



PYRESCOM

Mas des Tilleuls
66680 Canohès . France
aivia@pyres.com
www.aivia.pyres.com
Tél : +33 (0)4 68 68 39 68
Fax : +33 (0)4 68 68 39 69

